

FOR THE USE OF DENTAL PROFESSIONALS

# LS 360<sup>®</sup>

Lubricant Spray

## PRODUCT DESCRIPTION

**LS 360** is an oil spray used for cleaning and lubrication of all turbines, handpieces, contra-angles & air motors. It is specially formulated with paraffinic oils conforming to the highest standards of purity and safety as per pharmaceutical and medical device standards.

**LS 360** is designed to be used as an aerosol spray. It includes the actuator assembly that contains:

1. Spray Head
2. Spray Head Cap
3. Extension Nozzle
4. Pointed Adapter with Loop

**LS 360** is a unique and thermally stable material. Performing regular maintenance with **LS 360** improves the performance and extends the life time of the instruments.

Before handling, read product safety data sheet (SDS) and container labels for safe use.

## INTENDED USE / INDICATIONS FOR USE

For cleaning and lubrication of all turbines, handpieces, contra-angles & air motors.

## CONTRAINDICATIONS

Not for oral or topical use.

## DIRECTIONS FOR USE

- Lubricate prior to every sterilization.
- Hold the handpiece firmly wrapping it around a tissue paper.
- Hold the can in vertical position and spray into the handpiece for 2 to 3 seconds to clean and lubricate.
- Repeat until all debris is removed from the handpiece.
- Allow excess oil to drain by vertically positioning the handpiece head down.
- Sterilize after the oil is drained out completely.

## SALIENT FEATURES

- High purity & odorless.
- Optimized pressure with jet throw to ensure removal of debris with minimum material.
- Very easy to handle.
- Ensures longer life of instruments.
- Color stable and hydrophobic for longer shelf life.

## PRECAUTIONS

- Flammable. Pressurized container.
- Do not expose to temperatures exceeding 50°C.
- Do not pierce or burn even when apparently empty.
- Do not spray on a naked flame or any incandescent.
- Keep away from direct sunlight and heat sources.
- Do not spray near eyes or face.

- Avoid direct inhalation. Prolonged or frequent use of product may cause headache, drowsiness or skin irritation.
- Use in well ventilated area. Build-up of explosive mixture is possible without sufficient ventilation.
- Do not store in vehicle.

## STORAGE

The optimal storage condition is between 10°C to 30°C in a dry place.

Storage other than recommended temperature and conditions may affect the properties of the product.

## SHELF LIFE

When stored between 10°C and 30°C, this product has a shelf life of 36 months from the month of manufacturing.






## PACKING / CONTENT

The product is packed in 500 ml aerosol can and includes actuator assembly.







## NOTE

The material is for the use of dental professionals only. Processing should be carried out strictly according to the instructions for use. Liability cannot be accepted for damages resulting from misuse or failure to observe the instructions including storage conditions. The user is responsible for testing the product for its suitability and use for any other purpose that is not explicitly stated in the instructions.

## SYMBOLS

	Medical Device
	Lot No.
	Instruction for Use
	Storage Temperature
	Keep away from sunlight

Version / Date: 01 / November, 2024

	Manufactured by
	Date of Manufacture
	Date of Expiry
	Caution
	Flammable
	Ozone Friendly

**Anabond Stedman  
Pharma Research Pvt. Ltd.**  
No.117, Kattur Village, Vembudu Panchayat,  
Thirupurur Taluk, Chengalpattu District - 603 108,  
Tamil Nadu

ISO 13485:2016 Certified Company  
[www.anabondstedman.com](http://www.anabondstedman.com)  
For feedback / queries  
Email: [info@anabondstedman.com](mailto:info@anabondstedman.com)  
Mobile: +91 73059 65575

© - Registered Trademark

AS0002/01